

Zerlegbare Laminektomiezange Detachable laminectomy rongeurs


Durch Betätigung der einfachen, zum Patent angemeldeten Verriegelung, lässt sich der Schaft aufklappen und ermöglicht so die optimale Reinigung und Aufbereitung der Zangen.

Bei geöffnetem Schaft sind beide Teile noch über das Gelenk am Maulteil verbunden, wodurch verhindert wird, dass Einzelteile während des Aufbereitungsprozesses verwechselt werden oder verloren gehen.

Diese zukunftsweisende Technik wird den wachsenden Anforderungen im Aufbereitungsprozess gerecht und besticht durch intuitive und einfachste Handhabung.

By simply pressing the release button (patent pending) the shaft can be opened up, which enables the optimal cleaning and processing of the punches. Also in opened condition, both parts are still connected to each other in the jaw joint, which avoids loss and mistakes during the cleaning and sterilization process.

This forward-looking technology fulfills the growing demands related to the maintenance processes and impresses by its intuitive and simple handling.




AiM GmbH
Haldenstraße 9
78166 Donaueschingen
Germany

Telefon: +49 771 897807-0
Fax: +49 771 897807-33
info@aizm.de

Produktcharakteristik

Product characteristic

CUSHING



	2 mm x 10 mm	2 mm x 10 mm	2 mm x 10 mm
x = 130 mm - 5"	28030-13Z	28031-13Z	28033-13Z
x = 180 mm - 7"	28030-18Z	28031-18Z	28033-18Z

LOVE-GRUENWALD

	3 mm x 10 mm	3 mm x 10 mm	3 mm x 10 mm
x = 130 mm - 5"	28034-13Z	28035-13Z	28037-13Z
x = 180 mm - 7"	28034-18Z	28035-18Z	28037-18Z

SPURLING

	4 mm x 10 mm	4 mm x 10 mm	4 mm x 10 mm
x = 130 mm - 5"	28038-13Z	28039-13Z	28041-13Z
x = 180 mm - 7"	28038-18Z	28039-18Z	28041-18Z

AiM GmbH

Haldenstraße 9
78166 Donaueschingen
Germany

Telefon: +49 771 897807-0
Fax: +49 771 897807-33
info@aizm.de